

Biotech CE dialysis tubing

PRODUCT SPECIFICATIONS:

Biotech CE dialysis tubing	
Properties	Mildly hydrophobic, symmetric porosity
9 MWCO's	0.1-0.5 kD, 0.5-1.0 kD, 3.5-5 kD, 8-10 kD, 20 kD, 50 kD, 100 kD, 300 kD, 1000 kD
4 Flat Widths	100 mm, 16 mm, 24 mm, 31 mm
Quantity & Packaging	Wet in 0.05% sodium azide, 10 m/roll or 1 m/trial kit (includes 2 closures)

BIOTECH CE - Spectra Por® dialysis TUBING (10 meter/roll)

MWCO (kD)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	TUBING (10 m/roll) Part No.
0.1 - 0.5 kD	10	6.4	0.32	131048
	16	10	0.79	131054
	24	15	1.8	131057
	31	20	3.1	131060
0.5 - 1.0 kD	10	6.4	0.32	131084
	16	10	0.79	131090
	24	15	1.8	131093
	31	20	3.1	131096
3.5 - 5 kD	10	6.4	0.32	131192
	16	10	0.79	131198
	24	15	1.8	131201
	31	20	3.1	131204
8 - 10 kD	10	6.4	0.32	131264
	16	10	0.79	131270
	24	15	1.8	131273
	31	20	3.1	131276
20 kD	10	6.4	0.32	131336
	16	10	0.79	131342
	24	15	1.8	131345
	31	20	3.1	131348
50 kD	10	6.4	0.32	131372
	16	10	0.79	131378
	24	15	1.8	131381
	31	20	3.1	131384
100 kD	10	6.4	0.32	131408
	16	10	0.79	131414
	24	15	1.8	131417
	31	20	3.1	131420
300 kD	16	10	0.79	131450
	31	20	3.1	131456
1000 kD	16	10	0.79	131486
	31	20	3.1	131492

BIOTECH CE - Spectra Por® dialysis TRIAL KITS (1 meter/roll)

MWCO (kD)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	TUBING (1 m/roll) Part No.
0.1 - 0.5 kD	10	6.4	0.32	131048T
	16	10	0.79	131054T
	24	15	1.8	131057T
	31	20	3.1	131060T
0.5 - 1.0 kD	10	6.4	0.32	131084T
	16	10	0.79	131090T
	24	15	1.8	131093T
	31	20	3.1	131096T
3.5 - 5 kD	10	6.4	0.32	131192T
	16	10	0.79	131198T
	24	15	1.8	131201T
	31	20	3.1	131204T
8 - 10 kD	10	6.4	0.32	131264T
	16	10	0.79	131270T
	24	15	1.8	131273T
	31	20	3.1	131276T
20 kD	10	6.4	0.32	131336T
	16	10	0.79	131342T
	24	15	1.8	131345T
	31	20	3.1	131348T
50 kD	10	6.4	0.32	131372T
	16	10	0.79	131378T
	24	15	1.8	131381T
	31	20	3.1	131384T
100 kD	10	6.4	0.32	131408T
	16	10	0.79	131414T
	24	15	1.8	131417T
	31	20	3.1	131420T
300 kD	16	10	0.79	131450T
	31	20	3.1	131456T
1000 kD	16	10	0.79	131486T
	31	20	3.1	131492T