

## Biotech RC dialysis tubing

### PRODUCT SPECIFICATIONS:

<b>Biotech CE dialysis tubing</b>	
Properties	Mildly hydrophobic, symmetric porosity
9 MWCO's	0.1-0.5 kD, 0.5-1.0 kD, 3.5-5 kD, 8-10 kD, 20 kD, 50 kD, 100 kD, 300 kD, 1000 kD
4 Flat Widths	100 mm, 16 mm, 24 mm, 31 mm
Quantity & Packaging	Wet in 0.05% sodium azide, 10 m/roll or 1 m/trial kit (includes 2 closures)

### BIOTECH RC - Spectra Por® dialysis TUBING (5 meter/roll)

MWCO (kD)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	TUBING (5 m/roll) Part No.
0.1 - 0.5 kD	10	6.4	0.32	<a href="#">133192</a>
	16	10	0.79	<a href="#">133198</a>
8 - 10 kD	10	6.4	0.32	<a href="#">133264</a>
	16	10	0.79	<a href="#">133270</a>
20 kD	10	6.4	0.32	<a href="#">133336</a>
	16	10	0.79	<a href="#">133342</a>

### BIOTECH RC - Spectra Por® dialysis tRIAL KITS (0.5 meter/roll)

MWCO (kD)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	TUBING (0.5 m/roll) Part No.
3.5 - 5 kD	10	6.4	0.32	<a href="#">133192T</a>
	16	10	0.79	<a href="#">133198T</a>
8 - 10 kD	10	6.4	0.32	<a href="#">133264T</a>
	16	10	0.79	<a href="#">133270T</a>
20 kD	10	6.4	0.32	<a href="#">133336T</a>
	16	10	0.79	<a href="#">133342T</a>