

SpectraPor®

Float-A-Lyzer Dialysis Device

Membrane type	Biotech grade Cellulose Ester (CE)
Screw-on cap	Open and close for in-process sampling Color-coded for MWCO Polypropylene
Silicon O-ring	Leak proof and resealable
Flotation ring	Device body inserts through center to maintain buoyancy in buffer
Top/Bottom piece	Polyethylene
Potting	Polycarbonate
Sample loading	Polyurethane Disposable pipette for 5 and 10 ml (included) Pipette tip for 1 ml (not included)
Package/Quantity	Dry with glycerol, 12/package

9 MWCOS	COLOR CODE	QUANTITY
0.1 - 0.5 kD	Green	12/pkg
0.5 - 1 kD	Orange	12/pkg
3.5 - 5 kD	Black	12/pkg
8 - 10 kD	Yellow	12/pkg
20 kD	Red	12/pkg
50 kD	Violet	12/pkg
100 kD	Blue	12/pkg
300 kD	Amber	12/pkg
1000 kD	Pink	12/pkg

3 VOLUME SIZE	1 ml	5 ml	10 ml
Working volume	0.8 ml - 1.3 ml	4.0 ml - 5.5 ml	8.0 ml - 11.0 ml
Total length	5 cm	10 cm	16 cm
Membrane diameter	10 mm	10 mm	10 mm
Top piece diameter	23 mm	23 mm	23 mm
Flotation ring	38 mm	38 mm	38 mm