

SpectraPor®

Micro Float-A-Lyzer Dialysis Device

Membrane type	Biotech Grade Cellulose Ester (CE)
PP Luer-Lok cap	PP Luer-Lok cap Open and close for in-process sampling Color-coded for MWCO Leak proof and re-sealable
Body piece	Body piece Seals top/bottom ends of membrane Buoyant and self-orienting Polycarbonate
Potting	Potting Polyurethane
Sample loading	Sample loading 1 ml syringe (included)
Packaging	12/package

7 MWCOS	COLOR CODE	QUANTITY
0.1 - 0.5 kD	Green	12/pkg
0.5 - 1 kD	Orange	12/pkg
3.5 - 5 kD	Black	12/pkg
8 - 10 kD	Yellow	12/pkg
20 kD	Red	12/pkg
50 kD	White	12/pkg
100 kD	Blue	12/pkg

2 VOLUME SIZE	100-200 µl	400-500 µl
Height	4.4 cm	6.3 cm
Width	4.5 cm	4.5 cm
Thickness	1.9 cm	1.9 cm
Membrane flat width	10 mm	10 mm
Membrane diameter	6.4 mm	6.4 mm