

SpectraPor® (1-5)
DRY Standard Grade dialysis tubing

PRODUCT SPECIFICATIONS:

DRY Spectrapor® 1-5	
Properties	Mildly hydrophobic, symmetric porosity
9 MWCO's	0.1-0.5 kD, 0.5-1.0 kD, 3.5-5 kD, 8-10 kD, 20 kD, 50 kD, 100 kD, 300 kD, 1000 kD
4 Flat Widths	100 mm, 16 mm, 24 mm, 31 mm
Quantity & Packaging	Wet in 0.05% sodium azide, 10 m/roll or 1 m/trial kit (includes 2 closures)

Membrane Type	Purpose	MWCO	Flat Widths (mm)
SpectraPor® 1 Dialysis Membrane	Basic Dialysis	6 - 8 kD	10, 23, 32, 40, 50, 100, 120
SpectraPor® 2 Dialysis Membrane	Higher permeability	12 -14 kD	10, 25, 25, 105, 120
SpectraPor® 3 Dialysis Membrane	Basic Dialysis	3.5 kD	18, 45, 54
SpectraPor® 4 Dialysis Membrane	Basic Dialysis	12 -14 kD	10, 25, 32, 45, 75
SpectraPor® 5 Dialysis Membrane	Reinforced for weight	12 -14 kD	75, 100

SpectraPor® (1-5)
DRY Standard Grade dialysis tubing

DRY - Spectra Por® standard rc dialysis TUBING

View All	MWCO (kD)	Glycerin (yes/no)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	Length (m)	TUBING Part No.
Spectra Por® 1	6-8	yes	10	6.4	0.32	15	132645
		yes	23	14.6	1.7	30	132650
		yes	32	20.4	3.3	30	132655
		yes	40	25.5	5.1	30	132660
		yes	50	32	7.9	30	132665
		yes	100	64	32	15	132670
		yes	120	76	46	15	132675
Spectra Por® 2	12-14	yes	10	6.4	0.32	15	132676
		yes	25	16	2.0	15	132678
		yes	45	29	6.4	15	132680
		yes	105	67	34	15	132682
		yes	120	76	46	15	132684
Spectra Por® 3	3.5	yes	18	11.5	1.1	15	132720
		no	45	29	6.4	15	132724
		no	54	34	9.3	15	132725
Spectra Por® 4	12-14	yes	10	6.4	0.32	30	132697
		yes	25	16	2.0	30	132700
		yes	32	20.4	3.3	30	132703
		yes	45	29	6.4	30	132706
		yes	75	48	18	15	132709
Spectra Por® 5	12-14	no	75	48	18	15	132706
		noy	140	89	62	15	132709

DRY - Spectra Por® standard rc dialysis TRIAL KITS

View All	MWCO (kD)	Glycerin (yes/no)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	Length (m)	TRIAL KITS Part No.
Spectra Por® 1	6-8	yes	10	6.4	0.32	1,5	132645T1 , 132645T
		yes	23	14.6	1.7	1,5	132650T1 , 132650T
		yes	32	20.4	3.3	1,5	132655T1 , 132655T
		yes	40	25.5	5.1	1,5	132660T1 , 132660T
		yes	50	32	7.9	1,5	132665T1 , 132665T
		yes	100	64	32	1,5	132670T1 , 132670T
		yes	120	76	46	1,5	132675T1 , 132675T
Spectra Por® 2	12-14	yes	10	6.4	0.32	1,5	132676T1 , 132676T
		yes	25	16	2.0	1,5	132678T1 , 132678T
		yes	45	29	6.4	1,5	132680T1 , 132680T
		yes	105	67	34	1,5	132682T1 , 132682T
		yes	120	76	46	1,5	132684T1 , 132684T
Spectra Por® 3	3.5	yes	18	11.5	1.1	15	132720T1 , 132720T
		no	45	29	6.4	15	132724T1 , 132724T
		no	54	34	9.3	15	132725T1 , 132725T
Spectra Por® 4	12-14	yes	32	20.4	3.3	1,5	132703T1 , 132703T
		yes	75	48	18	1,5	132709T1 , 132709T
Spectra Por® 5	12-14	no	75	48	18	1,5	132754T1 , 132754T
		no	140	89	62	1,5	132757T1 , 132757T

Spectra Por® 1-4 READY-TO-USE DIALYSIS SACKS (for large volume)

View All	MWCO (kD)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	Volume (ml)	Qty.	1-4 READY-TO-USE Dialysis Sacks (for large volume)
Spectra Por® 1	6-8	50	32	7.9	80 - 400	10	132666
Spectra Por® 2	12-14	45	29	6.4	80 - 320	10	132681
Spectra Por® 3	3.4	54	34	9.3	90 - 450	10	132726
Spectra Por® 4	12-14	45	29	6.4	80 - 320	10	132707