

SpectraPor® (1-4)
DRY - Standard RC single-layer membrane

PRODUCT SPECIFICATIONS:

Membrane configuration	Dics	Square sheets
Properties	Hydrophilic, symmetric porosity	Hydrophilic, symmetric porosity
Size	33 mm, 47 mm & 100 mm dia	(refer to part number table)
Packaging	Dry with or without glycerin	Dry with or without glycerin
Quantity	50 discs/pkg	25 sheets/pkg

Membrane Type	Purpose	MWCO	Glycerin	Configurations
SpectraPor® 1	Basic Dialysis	6 - 8 kD	Yes	Discs and Sheets
SpectraPor® 2	Higher permeability	12 -14 kD	Yes	Discs and Sheets
SpectraPor® 3	Basic Dialysis	3.5 kD	No	Discs and Sheets
SpectraPor® 4	Basic Dialysis	12 -14 kD	Yes	Discs and Sheets

DRY - Spectra Por® 1-4 DIALYSIS MEMBRANE DISCS

View All	MWCO (kD)	Disc Diameter (mm)	Qty.	MEMBRANE DISC Part No.
Spectra Por® 1	6-8	33	50	132478
		47		132476
		100		132474
Spectra Por® 2	12-14	33	50	132482
		47		132480
		100		132477
Spectra Por® 3	3.4	33	50	132488
		47		132486
		100		132484
Spectra Por® 4	12-14	33	50	132498
		47		132496
		100		132494

DRY - Spectra Por® 1-4 DIALYSIS MEMBRANE DISCS

View All	MWCO (kD)	Disc Diameter (mm)	Qty.	MEMBRANE DISC Part No.
Spectra Por® 1	6-8	240 x 240	25	132677
Spectra Por® 2	12-14	200 x 200	25	132686
Spectra Por® 3	3.4	108 x 108	25	132723
Spectra Por® 4	12-14	150 x 150	25	132712