

SpectraPor® (6)

PREWETTED - Standard RC dialysis tubing

PRODUCT SPECIFICATIONS:

DRY Spectrapor® 1-5	
Properties	Mildly hydrophobic, symmetric porosity
9 MWCO's	0.1-0.5 kD, 0.5-1.0 kD, 3.5-5 kD, 8-10 kD, 20 kD, 50 kD, 100 kD, 300 kD, 1000 kD
4 Flat Widths	100 mm, 16 mm, 24 mm, 31 mm
Quantity & Packaging	Wet in 0.05% sodium azide, 10 m/roll or 1 m/trial kit (includes 2 closures)

PRE-WETTED - Spectra Por® 6 standard rc dialysis TUBING

MWCO (kD)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	TUBING Part No.
1 kD	18	11.5	1.1	132636
	38	24	4.6	132638
	45	29	6.4	132640
2 kD	18	11.5	1.1	132620
	38	24	4.6	132625
	45	29	6.4	132633
3.5 kD	18	11.5	1.1	132590
	45	19	6.4	132592
	54	34	9.3	132594
8 kD	12	7.5	0.45	132579
	18	11.5	1.1	128058
	24	15	1.8	132580
	32	20.4	3.3	132582
	40	25.5	5.1	132584
10 kD	12	7.5	0.45	132570
	18	11.5	1.1	128118
	24	15	1.8	132572
	32	20.4	3.3	132574
	45	29	6.4	132576
15 kD	12	7.5	0.45	132560
	18	11.5	1.1	128158
	24	15	1.8	132562
	32	20.4	3.3	132564
	45	29	6.4	132566
25 kD	12	7.5	0.45	132550
	18	11.5	1.1	128218
	28	18	2.5	132552
	34	22	3.7	132554
50 kD	10	6.4	0.32	132539
	12	7.5	0.45	132540
	28	18	2.5	132542
	34	22	3.7	132544

PRE-WETTED - Spectra Por® 6 standard rc dialysis TRIAL KITS

MWCO (kD)	Flat Width (mm)	Diameter (mm)	Vol/Length (ml/cm)	TRIAL KITS Part No.
1 kD	18	11.5	1.1	132636T
	38	24	4.6	132638T
	45	29	6.4	132640T
2 kD	18	11.5	1.1	132620T
	38	24	4.6	132625T
	45	29	6.4	132633T
3.5 kD	18	11.5	1.1	132590T
	45	19	6.4	132592T
	54	34	9.3	132594T
8 kD	12	7.5	0.45	132579T
	24	15	1.8	132580T
	32	20.4	3.3	132582T
	40	25.5	5.1	132584T
	50	32	7.9	132586T
10 kD	12	7.5	0.45	132570T
	24	15	1.8	132572T
	32	20.4	3.3	132574T
	45	29	6.4	132576T
	15 kD	12	7.5	0.45
24		15	1.8	132562T
32		20.4	3.3	132564T
45		29	6.4	132566T
25 kD		12	7.5	0.45
	28	18	2.5	132552T
	34	22	3.7	132554T
	50 kD	12	7.5	0.45
28		18	2.5	132542T
34		22	3.7	132544T